



TIMM23 Antibody

Product Code	CSB-PA023553GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14925
Immunogen	Human TIMM23
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	translocase of inner mitochondrial membrane 23 homolog (yeast);TIMM23;MGC22767;PRO1197;TIM23 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TIMM23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.