



TIMP1 Antibody

Product Code	CSB-PA023560GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01033
Immunogen	Human TIMP1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TIMP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.