



# TJAP1 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA023572GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q5JTD0   |
| <b>Immunogen</b>           | Human TJAP1  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, WB  |
| <b>Storage Buffer</b>      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.  |
| <b>Purification Method</b> | Antigen Affinity purified  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | tight junction associated protein 1<br>(peripheral); TJAP1; DKFZp686F06131; PILT; TJP4 ;   |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | TJAP1  |
| <b>Target Details</b>      | This gene encodes a tight junction-associated protein. Incorporation of the encoded protein into tight junctions occurs at a late stage of formation of the junctions. The encoded protein localizes to the Golgi and may function in vesicle trafficking. Alternatively spliced transcript variants have been described. A related pseudogene exists on the X chromosome. |