



TK1 Antibody

Product Code	CSB-PA023577GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04183
Immunogen	Human TK1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	thymidine kinase 1, soluble;TK1;TK2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.