



TNKS Antibody

Product Code	CSB-PA024004GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95271
Immunogen	Human TNKS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	tankyrase, TRF1-interacting ankyrin-related ADP-ribose polymerase;TNKS;PARP-5a;PARP5A;PARPL;TIN1;TINF1;TNKS1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNKS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.