



TOMM40 Antibody

Product Code	CSB-PA024050GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96008
Immunogen	Human TOMM40
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	translocase of outer mitochondrial membrane 40 homolog (yeast);TOMM40;C19orf1;D19S1177E;PER-EC1;PEREC1;TOM40 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TOMM40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.