



TPX2 Antibody

Product Code	CSB-PA024139GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW0
Immunogen	Human TPX2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	TPX2, microtubule-associated, homolog (Xenopus laevis);TPX2;C20orf1;C20orf2;DIL-2;DIL2;FLS353;GD:C20orf1;HCA519;HCTP4;REPP86;p100 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TPX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.