



TRAPPC4 Antibody

Product Code	CSB-PA024230GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y296
Immunogen	Human TRAPPC4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	trafficking protein particle complex 4;TRAPPC4;CGI-104;HSPC172;PTD009;SBDN;SYNBINDIN;TRS23 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRAPPC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.