



TRERF1 Antibody

Product Code	CSB-PA024411GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PN7
Immunogen	Human TRERF1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	transcriptional regulating factor 1; TRERF1; BCAR2; FLJ21198; HSA277276; RAPA; RP1-139D8.5; TReP-132; dJ139D8.5 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRERF1
Target Details	This gene encodes a zinc-finger transcriptional regulating protein which interacts with CBP/p300 to regulate the human gene CYP11A1.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.