



TRMT2B Antibody

Product Code	CSB-PA024540GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GJ1
Immunogen	Human TRMT2B
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	TRM2 tRNA methyltransferase 2 homolog B (S. cerevisiae); TRMT2B; CXorf34; FLJ12687; dJ341D10.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRMT2B
Target Details	This gene encodes a homolog of the TRM2 gene in <i>S. cerevisiae</i> . The yeast gene encodes a tRNA methyltransferase that plays a role in tRNA maturation. The yeast protein also has endo-exonuclease activity and may be involved in DNA double strand break repair. Alternative splicing results in multiple transcripts encoding different isoforms.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.