



TROAP Antibody

Product Code	CSB-PA025070GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12815
Immunogen	Human TROAP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	trophinin associated protein (tastin);TROAP;TASTIN ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TROAP