



TSHR Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA025131GA01HU |
| Storage | Store at -20°C. Stable for one year after shipment. Aliquoting is unnecessary for -20oC storage |
| Uniprot No. | P16473 |
| Immunogen | Human TSHR |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Tested Applications | ELISA,WB |
| Form | Liquid |
| Storage Buffer | PBS with 0.02% sodium azide and 50% glycerol pH 7.3 |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TSHR |