



TSR1 Antibody

Product Code	CSB-PA025183GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2NL82
Immunogen	Human TSR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	TSR1, 20S rRNA accumulation, homolog (S. cerevisiae);TSR1;FLJ10534;KIAA1401;MGC131829 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.