



TUBA8 Antibody

Product Code	CSB-PA025314GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY65
Immunogen	Human TUBA8
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.