



TUSC2 Antibody

Product Code	CSB-PA025351GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75896
Immunogen	Human TUSC2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	tumor suppressor candidate 2;TUSC2;C3orf11;FUS1;PAP;PDAP2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUSC2
Target Details	This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.