



TWIST1 Antibody

Product Code	CSB-PA025358GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15672
Immunogen	Human TWIST1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	twist homolog 1 (Drosophila);TWIST1;ACS3;BPES2;BPES3;SCS;TWIST;bHLHa38 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TWIST1
Target Details	Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. This protein is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.