



TXNDC12 Antibody

Product Code	CSB-PA025371GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95881
Immunogen	Human TXNDC12
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	thioredoxin domain containing 12 (endoplasmic reticulum);TXNDC12;AGR1;ERP18;ERP19;PDIA16;TLP19;hAG-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TXNDC12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.