



UBR2 Antibody

Product Code	CSB-PA025517GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWW8
Immunogen	Human UBR2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ubiquitin protein ligase E3 component n-recognin 2;UBR2;C6orf133;DKFZp686C08114;KIAA0349;MGC71112;RP3-392M17.3;bA49A4.1;dJ242G1.1;dJ392M17.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.