



UNC45B Antibody

Product Code	CSB-PA025626GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWX7
Immunogen	Human UNC45B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	unc-45 homolog B (C. elegans);UNC45B;CMYA4;FLJ38610;MGC119540;MGC119541;SMUNC45;UNC45 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UNC45B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.