



UNKL Antibody

Product Code	CSB-PA025647GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9P5
Immunogen	Human UNKL
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	unkempt homolog (Drosophila)-like;UNKL;C16orf28;FLJ12623;FLJ23360;KIAA0734;MGC5179;ZC3H5L;ZC3HD C5L ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UNKL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.