



USP46 Antibody

Product Code	CSB-PA025738GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62068
Immunogen	Human USP46
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ubiquitin specific peptidase 46;USP46;FLJ11850;FLJ12552;FLJ14283;FLJ39393 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	USP46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.