



UTP20 Antibody

Product Code	CSB-PA025766GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75691
Immunogen	Human UTP20
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	UTP20, small subunit (SSU) processome component, homolog (yeast);UTP20;DRIM ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UTP20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.