



VAC14 Antibody

Product Code	CSB-PA025779GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08AM6
Immunogen	Human VAC14
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Vac14 homolog (S. cerevisiae);VAC14;ArPIKfyve;FLJ10305;FLJ36622;FLJ46582;MGC149815;MGC149816;TAX1BP2;TRX ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VAC14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.