



VMA21 Antibody

Product Code	CSB-PA025866GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3ZQAQ7
Immunogen	Human VMA21
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	VMA21 vacuolar H ⁺ -ATPase homolog (S. cerevisiae);VMA21;MGC125514;MGC125516;MGC131652 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VMA21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.