



WDFY4 Antibody

Product Code	CSB-PA025993GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZS81
Immunogen	Human WDFY4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	WDFY family member 4;WDFY4;C10orf64;FLJ17416;FLJ36288;FLJ45748;KIAA1607;MGC40604;MGC45899 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	WDFY4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.