



YPEL1 Antibody

Product Code	CSB-PA026270GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60688
Immunogen	Human YPEL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	yippee-like 1 (Drosophila);YPEL1;FKSG3;MGC64992 ;
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	YPEL1
Target Details	This gene encodes a potential zinc-binding protein and is located in the region associated with DiGeorge syndrome on chromosome 22. The encoded protein may play a role in the regulation of cellular morphology and may function during craniofacial development.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.