



# YWHAH Antibody

<b>Product Code</b>	CSB-PA026289GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q04917
<b>Immunogen</b>	Human YWHAH
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, eta polypeptide;YWHAH;YWHA1 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	YWHAH
<b>Target Details</b>	This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and bovine orthologs. This gene contains a 7 bp repeat sequence in its 5 UTR, and changes in the number of this repeat have been associated with early-onset schizophrenia and psychotic bipolar disorder.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.