



ZBTB7A Antibody

Product Code	CSB-PA026343GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95365
Immunogen	Human ZBTB7A
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Protein A purified
Isotype	IgG
Alias	zinc finger and BTB domain containing 7A; ZBTB7A; DKFZp547O146; FBI-1; FBI1; LRF; MGC99631; ZBTB7; ZNF857A; pok emon ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZBTB7A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.