



ZC3H4 Antibody

Product Code	CSB-PA026363GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPT8
Immunogen	Human ZC3H4
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	zinc finger CCCH-type containing 4;ZC3H4;C19orf7;KIAA1064 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H4