



# ZEB1 Antibody

<b>Product Code</b>	CSB-PA026424GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P37275
<b>Immunogen</b>	Human ZEB1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	zinc finger E-box binding homeobox 1;ZEB1;AREB6;BZP;DELTA-EF1;MGC133261;NIL-2-A;NIL-2A;NIL2A;TCF8;ZEB;ZFHEP;ZFHX1A ;delta-EF1,
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ZEB1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.