



ZFYVE16 Antibody

Product Code	CSB-PA026470GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3T8
Immunogen	Human ZFYVE16
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	zinc finger, FYVE domain containing 16;ZFYVE16;DKFZp686E13162;ENDOFIN;KIAA0305 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZFYVE16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.