



ZKSCAN3 Antibody

Product Code	CSB-PA026494GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRR0
Immunogen	Human ZKSCAN3
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	zinc finger with KRAB and SCAN domains 3;ZKSCAN3;FLJ33906;KIAA0426;ZF47;ZFP306;ZNF306;ZNF309;ZSCAN13;Zfp47;dJ874C20.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZKSCAN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.