



ZNF395 Antibody

Product Code	CSB-PA026718GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8N7
Immunogen	Human ZNF395
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	zinc finger protein 395;ZNF395;DKFZp434K1210;HDBP2;PBF;PRF-1;PRF1;Si-1-8-14 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF395
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.