



ZXDC Antibody

Product Code	CSB-PA027161GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2QGD7
Immunogen	Human ZXDC
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ZXD family zinc finger C;ZXDC;DKFZp547N024;FLJ13861;FLJ36477;FLJ37830;MGC11349;ZXDL ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZXDC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.