



AFF1 Antibody

Product Code	CSB-PA001412GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51825
Immunogen	Human AFF1-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	AF4/FMR2 family, member 1; AFF1; AF-4; AF4; AF4-MLL; MGC134969; MLL/AF4; MLLT2; PBM1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AFF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.