



CHST13 Antibody

Product Code	CSB-PA005406GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NET6
Immunogen	Human CHST13-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	carbohydrate (chondroitin 4) sulfotransferase 13;CHST13;C4ST3;MGC119278;MGC119279;MGC119281 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHST13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.