



CLDN4 Antibody

Product Code	CSB-PA005506GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14493
Immunogen	Human CLDN4-specific
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	claudin 4;CLDN4;CPE-R;CPER;CPETR;CPETR1;WBSCR8;hCPE-R ;
Immunogen Species	Homo sapiens (Human)
Target Names	CLDN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.