



CPXM2 Antibody

Product Code	CSB-PA005930GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N436
Immunogen	Human CPXM2-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	carboxypeptidase X (M14 family), member 2;CPXM2;UNQ676 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CPXM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.