



MAPK1 Antibody

Product Code	CSB-PA013417GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28482
Immunogen	Human ERK1/2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	mitogen-activated protein kinase 1;MAPK1;ERK;ERK2;ERT1;MAPK2;P42MAPK;PRKM1;PRKM2;p38;p40;p41;p41mapk ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.