



# LGALS1 Antibody

<b>Product Code</b>	CSB-PA012882GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09382
<b>Immunogen</b>	Human GAL1, LGALS1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	lectin, galactoside-binding, soluble, 1; LGALS1; DKFZp686E23103; GAL1; GBP; Galectin-1, S-Lac lectin 1, HPL, HBL, HLBP14
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	LGALS1
<b>Target Details</b>	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.