



GPRASP1 Antibody

Product Code	CSB-PA009848GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5JY77
Immunogen	Human GPRASP1-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20?, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPRASP1