



BET1L Antibody

Product Code	CSB-PA002671GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYM9
Immunogen	Human GS15
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	blocked early in transport 1 homolog (S. cerevisiae)-like;BET1L;BET1L1;GOLIM3;GS15;HSPC197 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BET1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.