



HSD17B7 Antibody

Product Code	CSB-PA010776GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56937
Immunogen	Human HSD17B7P2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	hydroxysteroid (17-beta) dehydrogenase 7 pseudogene 2;HSD17B7P2;FLJ40773;HSD17B7;Hsd17b_2;MGC88684;bA291L22.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSD17B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.