



# HSP90AA1 Antibody

<b>Product Code</b>	CSB-PA010802GA01HU
<b>Abbreviation</b>	HSP90
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P07900
<b>Immunogen</b>	Human HSP90
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IF
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	heat shock protein 90kDa alpha (cytosolic), class A member 1;HSP90AA1;FLJ31884;HSP86;HSP89A;HSP90A;HSP90N;HSPC1;HSPCA;HSPCAL1;HSPCAL4;HSPN;Hsp89;Hsp90;LAP2 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	HSP90AA1