



# AIF1 Antibody

|                            |                                                                                                                                                                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA001490GA01HU                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                                     |
| <b>Uniprot No.</b>         | P55008                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>           | Human IBA1                                                                                                                                                                                                                                                                                                                        |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                            |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA,WB,IHC,IF                                                                                                                                                                                                                                                                                                                   |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                                                                                                                                                                                                                                              |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                                                                                                                                                                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                               |
| <b>Alias</b>               | allograft inflammatory factor 1;AIF1;AIF-1;IBA1;IRT-1 ;                                                                                                                                                                                                                                                                           |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                                                                                                                                                                                                                                                    |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                              |
| <b>Research Area</b>       | Neuroscience                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>        | AIF1                                                                                                                                                                                                                                                                                                                              |
| <b>Target Details</b>      | This gene is induced by cytokines and interferon. Its protein product is thought to be involved in negative regulation of growth of vascular smooth muscle cells, which contributes to the anti-inflammatory response to vessel wall trauma. Three transcript variants encoding different isoforms have been found for this gene. |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                                                                                                       |