



LRP2 Antibody

Product Code	CSB-PA013096GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98164
Immunogen	Human LRP2-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	low density lipoprotein-related protein 2;LRP2;DBS;gp330 ;Megalin
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.