



MUC16/CA125 Antibody

Product Code	CSB-PA015219GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI7
Immunogen	Human MUC16,CA125
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20?, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MUC16