



MYSM1 Antibody

Product Code	CSB-PA015382GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VVJ2
Immunogen	Human MYSM1-Specific
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Myb-like, SWIRM and MPN domains 1;MYSM1;2A-DUB;2ADUB;DKFZp779J1554;DKFZp779J1721;KIAA1915;RP4-592A1.1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MYSM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.