



MT-ND1 Antibody

Product Code	CSB-PA015076GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03886
Immunogen	Human ND1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NADH dehydrogenase, subunit 1 (complex I);ND1;MTND1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MT-ND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.