



NEO1 Antibody

Product Code	CSB-PA015712GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92859
Immunogen	Human NEO1-Specific
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	neogenin homolog 1 (chicken);NEO1;DKFZp547A066;DKFZp547B146;HsT17534;IGDCC2;NGN ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.